

AUTOCLAVE 1 - BI/0019 - CYCLE RESULTS

LEAK TEST REPORT:

CYCLE NUMBER: 60
USER NUMBER: 15
CYCLE SELECTED - P: P60
DATE: 21 11 2016
TIME: 11:00:47

| CYCLE PROGRESS: | TIME: | CHAMBER PRESS (kPa): | DRAIN TEMP ('C): |
|-------------------------|----------------|----------------------|------------------|
| Start: | 11:00:47 -78,0 | | 62,7 |
| Gasket Inflate: | 11:00:47 -78,0 | | 62,7 |
| Draw Vacuum: | 10:46:13 1,6 | | 62,4 |
| Settling: | 10:51:17 -80,0 | | 71,2 |
| Testing: | 10:53:17 -79,2 | | 66,8 |
| Air Admission: | 0,0 | | 0,0 |
| Gasket Release: | 0,0 | | 0,0 |
| COMPLETE: | | | |
| UNSTERILE: | 11:00:57 -78,0 | | 62,6 |
| STERILIZED FOR: min:sec | 0min0sec | | |
| CYCLE DURATION: min | 0min | | |

LOAD TEMP ('C): JACKET TEMP ('C):

| | |
|------|-------|
| 90,9 | 137,1 |
| 90,9 | 137,1 |
| 73,9 | 137,1 |
| 82,2 | 137,7 |
| 84,4 | 137,5 |
| 0,0 | 0,0 |
| 0,0 | 0,0 |